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Lumbar intervertebral disc replacement

The lumbar intervertebral disc replacement is made of technically pure titanium (Ti) or PEEK-OPTIMA™ HA Enhanced, with a osteoconductive surface. The osteoconductive surface of the implant ensures fast oseointegration (without any intermediary layer of fibrous tissue) and helps to create interbody fusion due to its osteoconductive properties. The osteoconductive intervertebral disc replacement is primarily designed to prevent the instability of a lumbar spine segment.

PLIF method

Characteristics: The lumbar implant has the basic shape of a prism or tapered prism and is made of technically pure titanium with a chemically modified osteoconductive surface (Fig. 1) or from PEEK-OPTIMA™ HA Enhanced by Invivo Ltd with a roughened osteoconductive surface (Fig. 2). The titanium implant is provided with two pairs of projections (wings) and the implant from PEEK-OPTIMA™ HA Enhanced has either ten pairs (21 mm long implant) or thirteen pairs (25mm long implant) of grooves based on the length of the implant. The grooves serve as stabilizing elements and are located on the sides of the implant. They adhere to the vertebral bodies after the implantation and contribute to the primary stability of the implant. The replacement made of PEEK-OPTIMA™ HA Enhanced is provided with three titanium markers for the visualization of the implant in X-ray examination. The IMPLASPIN implant is designed for lumbar intervertebral disc replacements to prevent instability of the affected spine segment (usually one or two segments). The implant is indicated for use only in skeletally mature patients.

Indication: Lumbar spine degenerative disease (primary and secondary instability, Grade I and II degenerative spondylolisthesis according to Meyerding)

Δ) *According to clinical findings (neurological examination, medical scoring systems, e.g. Oswestry Low Back Pain Disability Score, Visual Analogue Scale for Pain):*

Absolute indication: rapid progression of the paretic radicular lesion affecting both lower extremities, unilateral irritation-paretic radicular syndrome on the lower extremities, non-responding to conservative treatment, micturition disorder (cauda equina syndrome).

Relative indication: unilateral irritation radicular syndrome on lower extremities, non-responding to conservative treatment.

B) According to the imaging examination method:

MRI: content of the spinal canal and, especially, of the soft tissue structures, changes on intervertebral discs (acc. to Modic) and compression of the spinal nerve roots (intervertebral disc prolapse, cyst, hematoma, tumor, etc.).
CT: skeletal features of the affected spinal canal (joint osteoarthritis, lateral stenosis, etc.) and its content (compression of the spinal nerve roots by osteophytes).

Lumbar spine X-ray + dynamic imaging: degenerative changes, instability, skeletal spondylolisthesis grade I–II causing stenosis of the spinal canal, demonstration of intervertebral disc osteochondrosis (reduction of intervertebral disc height), demonstration of instability based on the development of anterior and posterior osteophytes and pathological movements during forward and backwards bending.

C) Electrophysiology:

EMG: identification of a radicular lesion with the development of a paretic component.

Contraindications: Contraindications may be relative or absolute. The surgeon must take these absolute and relative contraindications into account when making his/her decision. This list is not exhaustive.

Δ) *Absolute:* poor general health, patient's inability to undergo a longer surgical procedure due to serious reasons or severe disease, active systemic infection or infection localized to the site of implantation, neuromuscular disorder, allergy

or hypersensitivity to materials used: titanium and polyetheretherketone (PEEK).

B) *Relative:* any abnormality present which affects the normal process of bone remodeling (e.g. severe osteoporosis, poor bone mineral density), a specific/non-specific inflammatory lesion in the operated area, a tumorous osteolytic lesion in the operated area, metabolic disorder affecting osteogenesis, patient's inability to remember and comply with preoperative and postoperative instructions, conditions other than those indicated.

Adverse effects: *Possible adverse effects associated with general surgery* include, but are not limited to, complications from anesthesia, infection, hematoma, blood vessel damage, hemorrhage, nerve damage, organ or muscle damage, and death. *Possible adverse effects associated with lumbar spinal surgery* include, but are not limited to, neural structures damage (nerve root) with loss of sensation or mobility of the lower limbs and loss of bladder or bowel control, dural tear associated with CSF leakage, superficial or deep infection of the wound, postoperative instability, intrabdominal structures injury.

Possible adverse effects associated with the use of IMPLASPIN include: loosening or migration of the implant, postoperative instability, delayed union or non-union, subsidence of the implant into adjacent vertebral bodies, adjacent segment disease, loss of disc height, change of spine curvature, allergic reaction to the implant, failure of supplemental fixation. This list of adverse effects is not exhaustive. Additional complications that are not anticipated may also occur. Reoperation may be necessary to correct adverse effects.

Interactions: None observed to date.

Instructions for use: **A suitably chosen size** of the lumbar intervertebral replacement IMPLASPIN is inserted through the implant cut on the modified end of the insertion tool body and attached using the insertion screw passed through a hole in the insertion tool body and **manually** tightened (see the manual). The replacement from PEEK-OPTIMA™ HA Enhanced is inserted with a sliding motion to the desired position; the grooves are oriented into the vertebral bodies. The fastening screw is then released and the insertion tool body is removed. The titanium replacement is introduced in the direction the projections (wings) are oriented into the intervertebral space. After setting the titanium replacement into the desired position, a handle is connected to the end of the insertion tool body and it is stabilized against rotation using a fixation screw. Subsequently, the implant is rotated by 90° so that the projections (wings) point into the bodies of the adjacent vertebrae. Once the correct position is reached, the fixation screw is loosened and the handle is disconnected from the insertion tool body (see the manual). When choosing the correct size of the replacement, we recommend using IMPLASPIN tool – an implant analogue manufactured by LASAK.

Recommendations: *PLIF Surgery (Posterior Lumbar Interbody Fusion)* – IMPLASPIN is most often used for PLIF complemented with transpedicular fixation (Johnson & Johnson, Alinx, Novaspine, Stryker, etc.).

PLIF alone, without transpedicular fixation, can be used in the case of instabilities without translation shift.

Close contact with vital bone tissue and primary stability are required for undisturbed healing. On the second or third post-surgery day the patient is verticalized using prosthetic devices if needed (corset, crutches). The post-surgery clinical examination of the patient by a neurosurgeon and neurologist should be performed daily until complete healing of the surgical wound with additional regular follow-ups of the clinical condition 3–12, and 24 months afterwards. Spine X-ray and dynamic imaging should be performed concomitantly in two planes. This examination checks the fixation strength of the operated segment. Follow up CT/MRI examinations of the spinal canal should be performed at 3 and 9 months after the surgery. In these examinations, the extent of the decompression of the nerve structures in the operated segment and the contact of the implant with the bone tissue and possible signs of osteoblastic activity on the surface are determined.

TLIF method

Characteristics: The shape of the lumbar implant is based on the basic shape of a prism with a curved midline in the longitudinal direction. The implant is equipped with a gra-filed aperture connecting the sides of the prism that adhere to the

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vertebral bodies after the application, which contributes to the primary and secondary stability of the implant. All the edges of the cage are rounded in order to eliminate the risk of damage to nerve structures in the course of implantation. The implantation is always inserted unilaterally. The implant is made of technical-grade titanium with a chemically treated surface which ensures its osteoconductive properties. The osteoconductive surface of the implant ensures a firm bond with the bone tissue and exhibits osteoconductive properties. The color of the material is gray-black, its density is 4.500 kg/m³ and its tensile strength is at least 450 MPa. The implant is designed for intervertebral disc replacements to stabilize a lumbar spine segment affected by instability caused by a degenerative impairment of the spine (usually one or two segments) that must be fixed by a graft. The implant is indicated for use only in skeletally mature patients.

Indications: Degenerative disease of intervertebral discs (primary and secondary instability, Grade I and II degenerative spondylolisthesis according to Meyerding)

Δ) *According to a clinical finding (neurological examination, medical scoring systems, e.g. Oswestry Low Back Pain Disability Score, Visual Analogue Scale for Pain):*

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Relative indication: unilateral irritation radicular syndrome on lower extremities, non-responsive to conservative treatment.

B) According to imaging examination methods:

MRI: content of the spinal canal and, especially, of the soft tissue structures, changes on intervertebral discs (acc. to Modic) and compression of the spinal nerve roots (intervertebral disc prolapse, cyst, hematoma, tumor, etc.).

CT: skeletal features of the affected spinal canal (joint osteoarthritis, lateral stenosis, etc.) and its content (compression of the spinal nerve roots by osteophytes).

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Contraindications: May be relative or absolute. The surgeon must take these absolute and relative contraindications into account when making his/her decision. This list is not exhaustive.

Δ) *Absolute:* poor general health, patient's inability to undergo a longer surgical procedure due to serious reasons or severe disease, active systemic infection or infection localized to the site of implantation, neuromuscular disorder, allergy or hypersensitivity to materials used: titanium.

B) *Relative:* any abnormality present which affects the normal process of bone remodeling (e.g. severe osteoporosis, poor bone mineral density), a specific/non-specific inflammatory lesion in the operated area, a tumorous osteolytic lesion in the operated area, metabolic disorder affecting osteogenesis, patient's inability to remember and comply with preoperative and postoperative instructions and conditions other than those indicated.

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Possible adverse effects associated with the use of IMPLASPIN include: loosening or migration of the implant,

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