Instrumentarium IMPLASPIN®

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IMPLASPIN instrument set

Characteristics: The IMPLASPIN instrument set is a system designated for the application of IMPLASPIN intervertebral disc replacements. The instrument set is made of Ti6Al4V biocompatible titanium alloy according to ISO 5832-3 or of steel according to ISO 7153-1. The material has been certified for the specified application. The instrument set is intended for repeated use; however, it must be sterilized before each use.

Method of use:

Instrument set (ACIF) designated for IMPLASPIN cervical intervertebral disc replacement: The aperture of the implant is fitted on the end of the insertion tool body adjusted for the aperture. The implant is then fixed with the fastening screw which is pushed through the opening in the insertion tool body and firmly tightened by hand. After the insertion of the implant to the desired position, loosen the fastening screw and remove the insertion tool body.

Instrument set (PLIF) designated for IMPLASPIN lumbar intervertebral disc replacement: The aperture of the implant is fitted on the end of the insertion tool body adjusted for the aperture. The implant is then fixed with the fastening screw, which is pushed through the opening in the insertion tool body and firmly tightened by hand. After the insertion of the implant to the desired position, a handle is attached to the end of the insertion tool body and secured it against rotation using a fixation screw. The implant is then rotated about 90° so that the projections (wings) of the implant are directed to the adhering vertebral bodies. After reaching the position, the fixation screw is loosened and the handle is removed from the insertion tool body. Subsequently the fastening screw is loosened and the insertion tool body removed.

Instrument set (TLIF) designated for IMPLASPIN lumbar intervertebral disc replacement: The aperture of the implant is fitted on the end of the insertion tool body adjusted for the aperture. The implant is then fixed with the fastening screw, which is pushed through the opening in the insertion tool body and firmly tightened by hand. After the insertion of the implant the fastening screw is loosened and the insertion tool body is removed. Subsequently, the final insertion tool with one end adjusted for the shape of the implant, may be used for the desired positioning of the implant.

Please, contact the manufacturer in the event of any defects of the instrument set.

Warning: The instruments are designed exclusively for the application of IMPLASPIN implants – cervical and lumbar intervertebral disc replacement using a surgical technique provided in the instrument specification.

Packaging: The instrument set is supplied non-sterile. The instrument set is packaged in a STERIKING transparent pouch with indicators for steam sterilization and in a protective bag. The packages allow the instrument set to be removed from the pouch aseptically. The instrument set is intended for repeated use. The instrument set must not be used if the pack is damaged and must be disposed of. Disposal is not subject to special requirements.

Sterilization: Stainless steel and titanium surgical instruments are sterilized using moist heat in an autoclave (steam sterilizer). Other sterilization methods (dry heat, chemical) are not permitted. The sterilization of the instruments in an autoclave is performed in sterilization bags (paper/foil etc., use the types suitable for steam sterilization). The instruments are placed in the bags individually. The bags must not be piled in an unorganized fashion. To ensure proper steam penetration and adequate drying, use appropriate stands (racks). Prior to the first sterilization of new tools, it is necessary to remove the outer PE protective bag (sterilization is performed only in the primary – sterilization foil). Any other way of steam sterilization is not permissible (freely positioned unpackaged instruments, in mutual contact, etc.). Stainless steel instruments should be sterilized separately from the tools made of other materials. The maximum permitted load of the sterilizer must not be exceeded. Programs at 121 °C (exposure time of at least 20 min.) or 134 °C (exposure time of at least 4 min.) are suitable for the steam sterilization of stainless steel and titanium instruments. Sterilization temperatures higher than recommended can damage the instruments. Appropriate programs containing a drying phase must be used for the sterilization of packed instruments, (storage of instruments in wet sterilization bags can cause corrosion). Sterile instruments should be kept in the sterilization bags (protection against subsequent microbial contamination) at room temperature in a dry, dust-free and disinfected environment. Packages containing sterilized instruments must show the date of sterilization. If the sterilization period has expired, the instruments should be re-sterilized. It is necessary to monitor and record the sterilization cycle according to relevant regulations. Follow the instructions of the

sterilizer manufacturer and the applicable national and international legislation for instrument sterilization.

Cleaning and disinfection: Disinfection and cleaning of surgical instruments must be performed immediately after use. Organic residues on the instruments must not be left to dry on the instrument. Use only cleaning agents and disinfectants recommended for the surgical instruments. Disinfectants or cleaning agents containing a high percentage of chlorine, oxalic acid or hydrogen peroxide are not suitable for the cleaning and disinfection of stainless steel and titanium instruments. Disinfectants and cleaning products must be used in accordance with the manufacturer's instructions for use (as regards to concentration, reaction time and temperature). When disinfecting and cleaning the instruments use adequate protection.

The instrument surface is mechanically cleaned using a nylon brush (wire brushes may damage the surface of the instruments). Afterwards, an ultrasonic bath is used, where it is necessary to ensure that individual tools are not in contact. Particular attention should be paid to holes and cavities. Complex instruments should be cleaned in a dismantled state. The saline solutions, disinfection and cleaning agents residues must be rinsed thoroughly under running water and the instruments must be dried. Any contact between the instruments must also be avoided during the drying process. The drying temperatures should not exceed 120° C (temperatures too high may damage the instruments). Long exposure of the instrument to saline solution or evaporation of saline solution on the instrument's surface may cause its corrosion. Any contact between the instruments during cleaning or sterilization can cause corrosion of the instrument surface as well as any contact with an already corroded instrument. Cleaned and disinfected instruments must be kept in a dry place in the time between the surgeries (moisture can cause instrument corrosion). Damaged or corroded instruments must be discarded, further use is not permitted. Corroded instruments must not be sterilized (risk of contamination of the autoclave or other tools with rust particles and subsequent corrosion spread).

Storage: The instruments should be stored at room temperature and in a dry place.

Manufacturer: LASAK s.r.o., Českosobrodská 1047/46, 190 01 Prague 9 – Hloubětín, Czech Republic, www.lasak.com

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IMPLASPIN cervical gauges instrument set

Characteristics: IMPLASPIN cervical gauges is a system designated for the correct selection of IMPLASPIN cervical intervertebral disc replacements. IMPLASPIN cervical gauges are made of steel according to ISO 7153-1. The material has been certified for the specified application. The instrument set is intended for repeated use; however, it must be sterilized before each use.

Method of use: IMPLASPIN cervical gauges are inserted with a sliding movement in to the place where the future IMPLASPIN implant is to be positioned. The thickness of the gauge should be estimated according to the size of the intervertebral space to allow the appropriate thickness of the IMPLASPIN intervertebral implant to be determined. Please, contact the manufacturer in the cause of any defects of the cervical gauges.

Warning: IMPLASPIN cervical gauges are designated to be used exclusively for the selection of the optimum thickness of an IMPLASPIN cervical intervertebral disc replacement.

Packaging: The instrument set is supplied non-sterile. The instrument set is packaged in STERIKING transparent pouch with indicators for steam sterilization and in a protective bag. The packages allow the instrument set to be removed from the pouch aseptically. The instrument set is intended for repeated use. The instrument set must not be used if the pack is damaged and must be disposed of. Disposal is not subject to special requirements.

Sterilization: Stainless steel instruments are sterilized using moist heat in an autoclave (steam sterilizer). Other sterilization methods (dry heat, chemical) are not permitted. The sterilization of the instruments in an autoclave is performed in sterilization bags (paper/foil etc., use the types suitable for steam sterilization). The instruments are placed in the bags individually. The bags must not be piled in an unorganized

fashion. To ensure proper steam penetration and adequate drying, use appropriate stands (racks). Prior to the first sterilization of new tools, it is necessary to remove the outer PE protective bag (sterilization is performed only in the primary – sterilization bag). Any other way of steam sterilization is not permissible (freely positioned unpackaged instruments, in mutual contact, etc.). Stainless steel instruments should be sterilized separately from tools made of other materials. The maximum permitted load of the sterilizer must not be exceeded. Programs at 121 °C (exposure time of at least 20 min.) or 134 °C (exposure time of at least 4 min.) are suitable for the steam sterilization of stainless steel instruments. Sterilization temperatures higher than recommended can damage the instruments. Appropriate programs containing a drying phase must be used for the sterilization of packed instruments, (storage of instruments in wet sterilization bags can cause corrosion). Sterile instruments should be kept in the sterilization bags (protection against subsequent microbial contamination) at room temperature in a dry, dust-free and disinfected environment. Packages containing sterilized instruments must show the date of sterilization. If the sterilization period has expired, the instruments should be re-sterilized. It is necessary to monitor and record the sterilization cycle according to relevant regulations. Follow the instructions of the sterilizer manufacturer and the applicable national and international legislation for instrument sterilization.

Cleaning and disinfection: Disinfection and cleaning of surgical instruments must be performed immediately after use. Organic residues on the instruments must not be left to dry on the instrument. Use only cleaning agents and disinfectants recommended for the surgical instruments. Disinfectants or cleaning agents containing a high percentage of chlorine, oxalic acid or hydrogen peroxide are not suitable for the cleaning and disinfection of stainless steel instruments.

Disinfectants and cleaning products must be used in accordance with the manufacturer's instructions for use (as regards to concentration, reaction time and temperature). When disinfecting and cleaning the instruments use adequate protection.

The instrument surface is mechanically cleaned using a nylon brush (wire brushes may damage the surface of the instruments). Afterwards, an ultrasonic bath is used, where it is necessary to ensure that individual tools are not in contact. Particular attention should be paid to holes and cavities. Complex instruments should be cleaned in a dismantled state. The saline solutions, disinfection and cleaning agents' residues must be rinsed thoroughly under running water and the instruments must be dried. Any contact between the instruments must also be avoided during the drying process. The drying temperatures should not exceed 120 °C (temperatures too high may damage the instruments). Long exposure of the instrument to saline solution or evaporation of saline solution on the instrument's surface may cause its corrosion. Any contact between the instruments during cleaning or sterilization can cause corrosion of the instrument surface as well as any contact with an already corroded instrument. Cleaned and disinfected instruments must be kept in a dry place in the time between the surgeries (moisture can cause instrument corrosion). Damaged or corroded instruments must be discarded, further use is not permitted. Corroded instruments must not be sterilized (risk of contamination of the autoclave or other tools with rust particles and subsequent corrosion spread).

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